EN - Perineum Compress
DE - Perineumkompresse
FR - Compresses de périnée
NL- Perineum kompres
AR - جهاز ضغط العجان
CS - Komprese perinea
ES - Compresa de perineo
FI - Väililihaa kompressi
GR - Επίθεμα Περινέου
IT - Cuscinetto a Compressione per il Perineo
LT - Tarpvietės kompresas
NO - Perineum koprimerer
PL - Kompres na Krocze
PT - Compressa Para Períneo
RO - Compresă pentru perineu
RS - Perineum Komprese
RU - Компрессы для области промежности
SE - Perineum kompress
SK - Kompresia perinea

EN – Instructions for Compresses Warm I Cold

Cold Therapy - Freezer: Place the compress without cover (in a clean freezer bag) in the freezer for at least 2 hours before use.
Cold Water Bath: Submerge the compress in ice-cold water until the desired temperature is reached.

Warm Therapy

Do not leave the compress unattended during heating. Excessive heating can damage the compress.

Microwave: Heat the compress as indicated in the microwave. Then heat for 5 seconds each time to reach the desired temperature. Never heat for longer than 30 seconds.
Warm Water Bath: Bring water to a boil and remove the pan from the heat source. Place the compress in the warm water for 3-5 minutes (not boiling!).

Usage instructions: Dry the compress before use with a cloth and always check for damage or leaks before use. Insert the compresses into the provided covers when the desired temperature is reached. Note: always test the temperature before use. The compress can be used for a maximum of 20 minutes. After use, clean the compress with water and soap. Wash the covers at a maximum of 40 degrees Celsius.

Store the compress in a cool, dry place, away from sunlight.

Warning: Always follow these instructions for your safety and keep them. Keep out of reach of children and animals. CAUTION: Do not exceed the recommended heating time. Always remove the microwave rack and use the rotating glass plate; never place the compress on a metal surface. Always consult your microwave manual. Overheating can cause damage to the compress and surroundings. Check for damage before use; if damaged, discard immediately. Do not apply to open wounds or sensitive skin. Do not fall asleep while using the compress. Consult a doctor for extreme, prolonged, or chronic pain. For external use only, do not ingest the contents. Do not place the covers in the microwave. WonderMom is not liable for any damage resulting from incorrect use or failure to follow the instructions for use of the compresses.

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